

LabNotebook

Bijan Varjavand

November 3, 2016

1 Introduction

This day our group tested the hardness of all of our samples which we cut from the cold roller.

1.1 Goals

The goal was both to learn about using the Rockwell hardness tester as well as the properties that it measured.

2 Purpose

A crucial property of materials is their hardness, and learning about machinery and the value itself enhances the qualifications of each student as a materials scientist.

3 Setup

The samples for testing were recovered from storage the week before.

3.1 Materials

1018 Steel annealed, 6061 Al annealed.

3.2 Tools

A machine in the lab that heats and compresses polymers around a sample to form a puck. Rockwell hardness tester.

4 Procedure

We loaded each sample in a puck side by side, except for our thinnest sample which was loaded with another group's puck due to lack of space on our own.

After the pucks were made, their hardness was measured with the Rockwell machine (after calibration).

5 Results

Calibration data of a control piece of Al:

Sample Hardness(HRB)	57.0	57.8	57.6	57.0	58.3
----------------------	------	------	------	------	------

Hardness values can be found in the Lab Report.

6 Observations

Hardness of the materials changed with cold rolling - more cold rolled samples exhibited higher hardness across both Al and Steel.